

GentleCath Air™ for Women

URINARY CATHETER

Indications

GentleCath Air™ for Women intermittent catheter is a hydrophilic urinary catheter. It is intended for intermittent drainage of the urinary bladder of females who need assistance with drainage due to conditions causing urinary retention or dysfunction of the urinary system.



Duration of Use

The catheters have been designed for routine intermittent use. Devices are single-use only and must be removed immediately upon drainage of the bladder and disposed of. The overall usage of these devices can be expected to be ongoing repeated multiple times a day.

List of Product Range:

FR Size	Outer Diameter (mm)	Effective Length: 89.8mm (3.54")
		Overall Length: 134.8mm (5.307")
		REF
10	3.3	510600
12	4.0	510601
14	4.7	510602

Sterilization:

Sterile	Yes
Sterilized by	Irradiation
Shelf Life	24 months

Certification:

CE-Mark	Yes (EU 2797 UK 0086)
Product Classification	Class I Sterile (EU)
Quality System Certification	ISO 13485
Global Medical Device Nomenclature (GMDN)	45603

Material Specification:

Latex	No
DEHP	No
PVC	No
Catheter	Polyolefin-based Elastomer (POBE)
Funnel	High Density Polyethylene (HDPE) & Ethylene Vinyl Acetate (EVA)
Wetting Chamber	Polypropylene (PP) & Sterile Purified Water
Catheter Case	Polypropylene (PP)
Sealing Rings	Ethylene Propylene Diene Monomer (EPDM)

List of Product Range:

REF	FR Size	Shelf Box		Shipping Carton		
		Qty (pcs)	Dimensions (mm)	Qty (pcs)	Dimensions (mm)	Volume (m ³)
510600	10	30	144 x 158 x 57	480	610 x 184 x 233	0.0262
510601	12	30	144 x 158 x 57	480	610 x 184 x 233	0.0262
510602	14	30	144 x 158 x 57	480	610 x 184 x 233	0.0262

Product Component	What is it Made of?	Recyclable	Recycling Code
Carton and Instructions for Use	Paper/Cardboard	Yes	N/A
Catheter Case	Polypropylene (PP)	Yes	5
Catheter Tubing	Polyolefin-based Elastomer (POBE)	Medical waste. Once used, dispose of in normal household waste	N/A
Catheter Funnel	High Density Polyethylene (HDPE) & Ethylene Vinyl Acetate (EVA)	Medical waste. Once used, dispose of in normal household waste	N/A

Storage and Handling:




Keep dry



Do not use if package has been opened or damaged



Do not reuse



Manufacturer: **Convatec Limited**

Address: 200 Crossing Boulevard
Suite 101
Bridgewater, New Jersey 08807

Phone: 800.422.8811

Web: www.convatec.com



convatec

Continence Care

For further information, please contact your Convatec Continence Care representative or your local distributor of quality healthcare products.

To request samples to try, visit www.convatec.com.