

The most trusted line of foam dressings  
gives rise to a new *blue*.



Hydrofera Blue READY-*Border*.<sup>®</sup>



Restore the body's natural balance with *blue*.

Hydrofera**BLUE**<sup>®</sup>

# Hydrofera Blue READY-Border® provides a perfect solution for secure wound protection and atraumatic dressing removal.


## Clinical Benefits:

- Non-cytotoxic antibacterial protection
- Gentle silicone adhesive for atraumatic dressing changes
- Does not inhibit fibroblasts or growth factors

Secure hold, yet gentle removal

Waterproof and breathable cover film

Intimate contact of foam with the wound bed provides better absorption and bioburden control



**Hydrofera Blue® helps wounds heal faster at less cost<sup>1</sup>**


In a recent retrospective study examining wound treatment data of more than 30,000 patients over a two-year period, Hydrofera Blue® dressings combined with a comprehensive Advanced Wound Management Program showed faster healing times and lower treatment costs.

- Compared with other commonly used dressings on patients who were not on an Advanced Wound Management Program, Hydrofera Blue® delivered over a 50% reduction in healing time and a 75% reduction in nursing labor cost.

**All Hydrofera Blue® dressings are effective against a broad spectrum of microorganisms (gram-negative, gram-positive and yeast) commonly found in wounds.<sup>2</sup>**



# Hydrofera READY-Border® conforms to the most difficult-to-dress wounds.



## Perforated Silicone Adhesive Border

- Waterproof and breathable cover film
- Secure hold, yet gentle removal
- Intimate contact of foam to the wound bed



Secure fit for difficult-to-dress areas and mobility



Diabetic foot ulcers



Acute wounds



Gentle on skin tears



Leg ulcers




Atraumatic removal

### Hydrofera Blue CLASSIC<sup>®</sup>



#### Product No Standard Dressings

Product No	Product Size	Quantity/Box	HCPCS Codes*
HB2214	2" x 2" (5 cm x 5 cm)	10	A6209
HB4414	4" x 4" (10.2 cm x 10.2 cm)	10	A6209
HB4514	4" x 5" (10.2 cm x 12.7 cm)	10	A6210
HB6614	6" x 6" (15.2 cm x 15.2 cm)	10	A6210
HBT0906	9-mm Tunneling Dressing (1.2 g)	10 pair	A6215



#### Heavy Drainage Dressings

HBHD4450	4" x 4" x 0.5" thick (10.2 cm x 10.2 cm x 1.2 cm)	5	A6209
HBHD6675	6" x 6" x 0.75" thick (15.2 cm x 15.2 cm x 1.9 cm)	5	A6210



#### Ostomy Dressing

HBRF2650	2.5" (6.4 cm) Diameter Ostomy Dressing	10	A6209
----------	--	----	-------

### Hydrofera Blue COMFORTCEL<sup>®</sup>



HBCC4550	4" x 5" x 0.5" (10.2 cm x 12.7 cm x 1.3 cm)	5	A6210
HBCC6675	6" x 6" x 0.75" (15.2 cm x 15.2 cm x 1.9 cm)	5	A6210

### Hydrofera Blue READY<sup>®</sup>



HBRS2520	2.5" x 2.5" (6.4 cm x 6.4 cm)	10	A6209
HBRS4520	4" x 5" (10.2 cm x 12.7 cm)	10	A6210
HBRS8820	8" x 8" (20.3 cm x 20.3 cm)	10	A6211



#### Tube Site Dressing (TSD)

HBR3535	3.5" x 3.5" (8.9 cm x 8.9 cm)	30	A6209
---------	-------------------------------	----	-------

### Hydrofera Blue READY-Border<sup>®</sup>



HBRB4040	4" x 4" (10.2 cm x 10.2 cm)	10	A6212
HBRB5959	5.9" x 5.9" (15 cm x 15 cm)	10	A6212
HBRB6080	6" x 8" (15 cm X 20.3 cm)	10	A6213

### Hydrofera Blue TRANSFER<sup>®</sup>



HBRT2525	2.5" x 2.5" (6.4 cm x 6.4 cm)	10	A6209
HBRT4050	4" x 5" (10.2 cm x 12.7 cm)	10	A6210
HBRT8080	8" x 8" (20.3 cm x 20.3 cm)	10	A6211

1. Hurd T. Improving the quality of chronic wound care using an advanced wound management program and Gentian Violet / Methylene Blue-Impregnated Antibacterial Dressings: A Retrospective Study. Surgical Technology International Volume 35, Sept 2019  
2. Data on file at Hydrofera, LLC.

The BSN Medical / Essity Inc. Reimbursement Disclaimer:

The Centers for Medicare & Medicaid Services (CMS) uses Palmetto GBA LLC., as the Pricing, Data Analysis, and Coding (PDAC) Contractor. The PDAC Contractor maintains the Product Classification List, located at <https://www4.palmettogba.com/pdac.dmecls/>. The Product Classification List is a searchable database containing products that have received HCPCS coding verification from PDAC and their corresponding HCPCS codes. BSN Medical Inc., however, makes no representations as to the accuracy of the information contained within the PDAC database, nor any representations as to whether its products are reimbursable under any government sponsored healthcare program and/or private-insurance program. It is the provider's sole responsibility, in consultation with the insurer, if necessary, to determine medical necessity, ensure coverage criteria is met, and submit appropriate codes, modifiers, and charges for services rendered. RV: 12/24

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician (or properly licensed practitioner). Prior to the use of this device be sure to read the complete information on how to use this device, including Warnings, Cautions and Instructions for Use and all other package inserts supplied with the product.