

## Navina Classic Product Information

### Product description

Transanal Irrigation System with manual control unit.

### Shelf life and uses

All Navina Irrigation System components have a shelf life of 3.5 years, except the Navina Catheters and the Navina Cone, which have 2 years. When Navina components are packed in sets, the component with shortest shelf life determine the shelf life of the set.

Navina catheters and cones are single use, the water container has 15 uses, and the Navina Classic control unit and tube set has 100 uses.

### Material composition

**Catheter Material:** Styrene Ethylene Butylene Styrene (SEBS)/ Polypropylene (PP) with a hydrophilic Polyvinylpyrrolidone (PVP) coating.

**Primary Package Material:** PE/PET Pouch.

**Secondary Package Material:** Cardboard and corrugated cardboard.

### Intended use

Navina Irrigation System is intended for Transanal Irrigation by instilling water up into the lower part of the colon through a rectal catheter.

### Indications for use

The Navina Irrigation System is indicated to help adults and children from 2 years who suffer from fecal incontinence, chronic constipation, and/or time consuming bowel management. By instilling water up into the lower part of the colon, the peristaltic muscles in the bowel can be triggered and start to evacuate the content of the lower colon and rectum.

### Contraindications

**Do NOT** use Navina Systems if you have one or more of the following:

- Known anal or colorectal stenosis.
- Active inflammatory bowel disease.
- Acute diverticulitis.
- Colorectal cancer.
- Ischemic colitis.
- You are within three months of anal or colorectal surgery.
- You are within 4 weeks of previous endoscopic polypectomy.

As the list may not be exhaustive, healthcare professionals will always consider individual user factors as well.



### Navina Classic System

Reference	Type	Specification
6910540	Navina Classic System	1 Navina Classic control unit 1 Water container 1 Tube set 1 Navina case 1 Accessory set
6910640	Navina Classic System Refill	1 Navina Classic control unit 1 Water container 1 Tube set
6900340 6900440 6901640	Navina Consumable Set Regular Navina Consumable Set Small Navina Consumable Set Cone	1 Water container 15 Rectal catheters
6894040 6894340 6894640	Navina Catheter Set Regular Navina Catheter Set Small Navina Catheter Set Cone	10 Rectal catheters
6901140	Navina Tube Set	1 Catheter tube 1 Water container tube
6901200	Navina Accessories	2 Grip rings (1 vertical, 1 horizontal) 1 Lanyard 1 Positioning strap 1 Navina Smart positioning clip
6901340	Navina Classic control unit	1 Navina Classic control unit

**wellspect.us**

Wellspect HealthCare, 1235 Friendship Rd, Suite 205, Braselton, GA 30517.  
Toll-free: +1 877 456 3742 .

**Navina™ Classic**



Manufacturer: Wellspect HealthCare, Aminogatan 1, P.O. Box 14, SE-431 21 Mölndal, Sweden.  
Phone: +46 31 376 40 00. [www.wellspect.us](http://www.wellspect.us)



Rx Only